Urgent Warning

Urgent device correction to the Continuous Renal Replacement Therapy device - Prismaflex issued by Baxter Healthcare.

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Saudi Food and Drug Authority SFDA would like to bring to your attention that Baxter Healthcare is initiating a field action in order to update software versions of the Prismaflex Control Unit - Continuous Renal Replacement Therapy. Baxter has received reports of device operators failing to adhere to the instructions for use pertaining to the safe unloading of disposable sets from the Prismaflex Control Unit. These steps are required to safely disconnect the patient before proceeding to unload the filter set after treatment. Unloading of the disposable set without following the instructions and warnings on the Prismaflex Control Units may lead to severe blood loss and potentially fatal outcomes.

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For further information, please see the recall by [Click Here](#).

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Recommendation:

- Contact Baxter AR to arrange for software update
- Prior to this software update and to ensure that the operators correctly unload the disposable set, Baxter is reminding its customers to follow the instructions provided in the Operator's manual, and displayed on screen, which are as follows:
 - 1- Clamp all lines,
 - 2- Disconnect Access and Return blood line from the blood access device, and
 - 3- Verify that all lines are clamped and the patient is disconnected.

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You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

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Please **share** this information with Healthcare Professionals in your facility.

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Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

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National Center for Medical Devices Reporting (NCMDR Website).

Medical Devices Sector
Saudi Food and Drug Authority
Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292)
North Ring Road - Al Nafal Unit (1)
Riyadh 13312 - 6288
Tel: +966 (11) 2038222 Ext: 2406 -2412 - 2430
Fax: +966 (11) 2757245

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For latest published Recalls/Alerts, please visit NCMDR website:

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