Safety Communication

Risk of occlusion of the airway when using Solus airway tube manufactured by Intersurgical Ltd.

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Problem:

Due to a manufacturing fault, inflation of the cuff may reduce the airway tube cross-section, leading to partial or total occlusion of the airway and reduced air flow to the patient.

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Recommendation:

- Identify affected devices (see below) and stop using them.
- Quarantine all affected devices in stock.
- Contact the manufacturer's authorized representative in KSA for corrective action.
- Share this information with Healthcare Professionals in your facility.

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The following models and LOT numbers are affected:

- Solus Flexible, small paediatric, wire-reinforced laryngeal mask airway 2
- Solus Flexible, large paediatric, wire-reinforced laryngeal mask airway 2.5
- Solus Flexible, small adult, wire-reinforced laryngeal mask airway 3
- Solus Flexible, medium adult, wire-reinforced laryngeal mask airway 4
- Solus Flexible, large adult, wire-reinforced laryngeal mask airway 5

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You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

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Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

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National Center for Medical Devices Reporting ([NCMDR Website](#)).

Medical Devices Sector

Saudi Food and Drug Authority

Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292)

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For latest published Recalls/Alerts, please visit NCMDR website:

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