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Malfunction in Intra-Aortic Balloon Pumps (IABP) manufactured by Maquet

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Problem:

A complaint was received involving a CS300 IABP that did not pump due to an electrical test failure code #58 (power up vent tests fail), maintenance code #3, and an autofill failure, which has been associated to a patient death due to the failure of the device to initiate therapy.

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Recommendation/Actions:

- Pursuant to the WARNINGS section of the Operating/User Instructions, clinicians are instructed not to leave the patient unattended during IABP therapy;
- Until the service is performed, Getinge recommends powering on the IABP prior to inserting the IAB catheter to allow the IABP to successfully complete its self-test. This action will

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
- take less than 60 seconds to perform;
- In the event the IABP fails to complete the self-test and exhibits electrical test failure code 58, remove the IABP from service and contact your local Maquet/Getinge Sales & Service Office;
 - In the unlikely event that a sudden interruption of therapy occurs, transfer the patient to an alternative IABP. If an alternative IABP is unavailable; manually inflate the IAB with air or helium and immediately aspirate. Please refer to the intra-aortic balloon catheter instructions for use; and
 - Each affected facility will be contacted by Maquet/Getinge Service Team to schedule on-site service.

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For further information, please see the Safety notice ([Click Here](#)).

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You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

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Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

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National Center for Medical Devices Reporting ([NCMDR Website](#)).

Medical Devices Sector

Saudi Food and Drug Authority

Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292)

North Ring Road - Al Nafal Unit (1)

Riyadh 13312 - 6288

Tel: +966 (11) 2038222 Ext: 2406

Fax: +966 (11) 2757245

For latest published Recalls/Alerts, please visit NCMDR website:

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