

## Safety Communication

### accelerated battery depletion cause elevated likelihood for early replacement

<b>Device/ Product Description:</b>	Subcutaneous Implantable Cardioverter Defibrillators
<b>Brand Name:</b>	EMBLEM MRI S-ICDs
<b>Lot numbers/Serials:</b>	213750, 213773, 213831, 213852, 213925, 213954, 213960, 213971, 214034, 214037, 214051, 214072
<b>Manufacturer:</b>	Boston Scientific
<b>Problem:</b>	performance of EMBLEM™ Subcutaneous Implantable Cardioverter Defibrillators (S-ICDs) may result in a need for device replacement (ERI/EOL) earlier than expected due to compromised performance of an electrical component causing accelerated battery depletion.
<b>Recommendation/Actions:</b>	<ul style="list-style-type: none"> <li>- monitor patients via LATITUDE to facilitate prompt detection of ERI/EOL</li> <li>- Perform a device follow-up every 3 months via remote or in-office interrogation.</li> <li>- Promptly investigate any suspected indication of accelerated depletion</li> <li>- Evaluate Risk for patients : <ul style="list-style-type: none"> <li>• with a history of life-threatening ventricular arrhythmias such as a secondary prevention</li> <li>• indication or previous appropriate shock for VT/VF2 .</li> <li>• ☒who are unable to be reliably followed every 3 months (via LATITUDE and/or in-clinic interrogation).</li> <li>• ☒who are not monitored via LATITUDE and are unable to hear beeping tones.</li> </ul> </li> <li>- Replace As Needed</li> </ul> <p>You can find more information and recommendations from ( <a href="#">HERE</a>).</p>



You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative of your product for corrective action.

Healthcare Professionals should **report any adverse events** suspected to be associated with affected devices above (or other Medical Devices) to:

**[National Center for Medical Devices Reporting.](#)**

Medical Devices Sector

Saudi Food and Drug Authority

Postal Address: Saudi Arabia - Saudi Food and Drug Authority

4904 northern ring branch rd - Hitteen Dist.

RIYADH 13513 - 7148

Tel: +966 (11) 2038222 Ext: 2995, 2952

Fax: +966 (11) 2757245

Or

**[Saudi Vigilance](#)**

For latest published Recalls/Alerts, please visit ([NCMDR Website](#))

Sincerely,  
NCMDR Team