

Safety Communication

رسالة سلامة

Issues related to the polyester-based polyurethane (PE-PUR) sound abatement foam used in Philips Continuous and Non- Continuous Ventilators

Device/ Product Description:	Continuous Ventilator	
Affected product:	All Devices manufactured before 26 April 2021, All serial numbers	
	Continuous Ventilator	Trilogy 100
		Trilogy 200
		Garbin Plus, Aeris, LifeVent
	Continuous Ventilator, Minimum Ventilatory Support, Facility Use	A-Series BiPAP Hybrid A30
		A-Series BiPAP V30 Auto
	Continuous Ventilator, Non-life Supporting	A-Series BiPAP A40
A-Series BiPAP A30		
Manufacturer:	Philips Respironics	
Problem:	<p>Issues can result in serious injury which can be life threatening, cause permanent impairment, and/or require medical intervention to preclude permanent impairment:</p> <ol style="list-style-type: none"> 1) polyester-based polyurethane (PE-PUR) foam may degrade into particles which may enter the device's air pathway and be ingested or inhaled by the user 2) The PE-PUR foam may off-gas certain chemicals. The foam degradation may be exacerbated by use of unapproved cleaning methods, such as ozone, and off-gassing may occur during operation. 	

<p>Recommendation /Actions:</p>	<ul style="list-style-type: none"> • Make sure that this document is reached to the end-users. • Ensure that all devices affected are identified. • Follow the instructions as indicated in the attached Field Safety Corrective Action. • Do not stop or alter a prescribed therapy before a discussion with the physician. At the discretion of the treating clinical team, the benefit of continued usage of these ventilator devices may outweigh the risks. • Use an inline bacterial filter If the physician determines to continue using this device, Consult the Instructions for Use for guidance on installation. • You can refer to below link for questions and answers: https://www.philips.sa/en/healthcare/e/sleep/communications/src-update • Contact the authorized representative for required corrective action. <p>For more information please click here.</p> <p>If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through: NCMDR Vigilance system (19999)unified call center</p>	
<p>Devices/Products photo:</p>	<p>Devices/Products Photos are available at the following link: https://www.philips.sa/en/healthcare/e/sleep/communications/src-update</p>	
<p>Authorized Representative Details</p>	<p>AR name:</p>	<p>Philips Healthcare Saudi Arabia Ltd.</p>
	<p>Assigned Contact Person:</p>	<p>Mohammed Alsamhan</p>
	<p>Mobile/Phone:</p>	<p>+966 503879145</p>
	<p>Email:</p>	<p>SFDA.SA.MET@Philips.com KSA.CCC@Philips.com</p>