

**Safety Communication**

**رسالة سلامة**

**Cardio Products exploitable remotely**

<b>Device/ Product Description:</b>	Q-Stress, X-Scribe, H-Scribe, R-Scribe	
<b>Affected product:</b>	H-Scribe	V 5.01 – V 6.4.0
	Q-Stress	V 6.0.0 – V 6.3.1
	X-Scribe	V 5.01 – V 6.3.1
	R-Scribe	V 5.01 – V 7.0.0
<b>Manufacturer:</b>	Hillrom - WelchAllyn Inc.	
<b>Problem:</b>	Authentication Bypass Using an Alternate Path or Channel, Successful exploitation of this vulnerability could allow an attacker to access privileged accounts.	

<p><b>Recommendation /Actions:</b></p>	<p>Hillrom plans to release software updates to address this vulnerability in their next software release.</p> <p>In the interim:</p> <ul style="list-style-type: none"> <li>- Hillrom recommends the following workaround and mitigation to reduce the risk:             <ul style="list-style-type: none"> <li>• Disable the SSO feature in the respective Modality Manager Configuration settings. Please refer to the instructions for use (IFU) and/or service manual for instructions on how to disable SSO.</li> </ul> </li>   <li>- Hillrom recommends users upgrade to the latest product versions when updated products are available. Information on how to update these products to their new versions can be found on the <a href="#">Hillrom disclosure page</a>.</li>   <li>- Hillrom recommends the following additional workarounds to help reduce risk:             <ul style="list-style-type: none"> <li>• Apply proper network and physical security controls.</li> <li>• Apply authentication for server access.</li> </ul> </li> </ul> <p>For more information, please click <a href="#">here</a>.</p> <p>If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through:  <a href="#">NCMDR</a>  <a href="#">Vigilance system</a>            (19999) unified call center</p>	
<p><b>Authorized Representative Details</b></p>	<p>AR name:</p>	<p>Farouk Mamoun Tamer &amp; Co.</p>
	<p>Assigned Contact Person:</p>	<p>Hamad Alnami</p>
	<p>Mobile/Phone:</p>	<p>+(966) 558121998</p>
	<p>Email:</p>	<p><a href="mailto:RAPMS@tamergroup.com">RAPMS@tamergroup.com</a></p>