

Medical Devices Sector  
National Center for Medical Devices Reporting

رسالة سلامة  
**Safety Communication**

To: Healthcare providers	إلى: مقدمي الرعاية الصحية
--------------------------	---------------------------

<b>Title:</b>	MiniMed™ 600 and 700 series insulin pump Battery Cap	العنوان:												
<b>Medical Device Description</b>	<table border="1"> <tr> <td>MiniMed™ 640G Insulin Pump</td> <td>MMT-1711, MMT-1712, MMT-1751, MMT-1752</td> </tr> <tr> <td>MiniMed™ 670G Insulin Pump</td> <td>MMT-1761, MMT-1762, MMT-1781, MMT-1782</td> </tr> <tr> <td>MiniMed™ 720G Insulin Pump</td> <td>MMT-1809, MMT-1810, MMT-1859, MMT-1860</td> </tr> <tr> <td>MiniMed™ 740G Insulin Pump</td> <td>MMT-1811, MMT-1812, MMT-1861, MMT-1862</td> </tr> <tr> <td>MiniMed™ 770G Insulin Pump</td> <td>MMT-1881, MMT-1882, MMT-1891, MMT-1892</td> </tr> <tr> <td>MiniMed™ 780G Insulin Pump</td> <td>MMT, 1885, MMT-1886, MMT-1895, MMT-1896</td> </tr> </table>	MiniMed™ 640G Insulin Pump	MMT-1711, MMT-1712, MMT-1751, MMT-1752	MiniMed™ 670G Insulin Pump	MMT-1761, MMT-1762, MMT-1781, MMT-1782	MiniMed™ 720G Insulin Pump	MMT-1809, MMT-1810, MMT-1859, MMT-1860	MiniMed™ 740G Insulin Pump	MMT-1811, MMT-1812, MMT-1861, MMT-1862	MiniMed™ 770G Insulin Pump	MMT-1881, MMT-1882, MMT-1891, MMT-1892	MiniMed™ 780G Insulin Pump	MMT, 1885, MMT-1886, MMT-1895, MMT-1896	اسم ووصف الجهاز/المستلزم الطبي:
MiniMed™ 640G Insulin Pump	MMT-1711, MMT-1712, MMT-1751, MMT-1752													
MiniMed™ 670G Insulin Pump	MMT-1761, MMT-1762, MMT-1781, MMT-1782													
MiniMed™ 720G Insulin Pump	MMT-1809, MMT-1810, MMT-1859, MMT-1860													
MiniMed™ 740G Insulin Pump	MMT-1811, MMT-1812, MMT-1861, MMT-1862													
MiniMed™ 770G Insulin Pump	MMT-1881, MMT-1882, MMT-1891, MMT-1892													
MiniMed™ 780G Insulin Pump	MMT, 1885, MMT-1886, MMT-1895, MMT-1896													
<b>Manufacturer:</b>	Medtronic	اسم المصنع:												
<b>Authorized Representative</b>	Medtronic Saudi Arabia	الممثل المعتمد												

<b>Medical Devices Marketing Authorization (MDMA)</b>	ME0000000248SFDAA00008, ME0000000248SFDAA00018, ME0000015320SFDAA 00001	إذن التسويق:
<b>Potential /Associated Risks</b>	A potential issue relating to the battery cap on the pumps. If the metal contact on the battery cap becomes loose or falls off from the battery cap, it can result in an incomplete battery connection, leading to no power source to the pump. When the pump detects no power source, an “Insert battery” alarm will occur, and insulin delivery will immediately stop. After 10 minutes, the alarm sound may increase to a siren, and the pump will turn off.	المخاطر المحتملة المرتبطة بالجهاز أو المستلزم الطبي
<b>Recommendations</b>	<p>Actions required by Physician, Healthcare Professional:</p> <ol style="list-style-type: none"> <li>1. Inform impacted pump users of the MiniMed™ 600 and/or 700 series insulin pump about the Field Safety Notice using the enclosed pump user letter.</li> </ol>	التوصيات

	<p>2. Please assist pump users in locating and inspecting the battery cap on their insulin pumps per the instructions provided in the pump user letter.</p> <p>Refer to the attached FSN for more details.</p> <p><a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=16159">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=16159</a></p>	
--	---	--

<b>For Reporting:</b>	  	للإبلاغ:
-----------------------	---	----------