Date: 10 Mar 2024

Reference Number: SG-2403-416-H

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Medical Devices Sector
National Center for Medical Devices Reporting

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To: Healthcare Provider		إلى: مقدمي الرعاية الصحية
Title	Recommendations for healthcare providers regarding the reprocessing procedures of duodenoscopes	العنوان
Medical Device Description	Reusable Duodenoscopes	اسم ووصف الجهاز/المستلزم الطبي
Manufacturer	All	اسم المصنع
Authorized Representative	All	الممثل المعتمد
Medical Devices Marketing Authorization (MDMA)	All	إذن التسويق
Potential /Associated risks	Based on the SFDA's post-market clinical evaluation study to analyze and evaluate the safety of reusable duodenoscope devices, the study findings indicated a potential risk of increasing contamination rate associated with reprocessed duodenoscopes.	المخاطر المحتملة/ المرتبطة بالجهاز أو المستلزم الطبي
Recommendations	<p>Recommendations for healthcare providers:</p> <ul style="list-style-type: none"> -Follow the manufacturer's instructions for reprocessing procedures of duodenoscopes. -Use appropriate cleaning equipment for duodenoscopes in compliance with the manufacturer's instructions for use. -Healthcare providers should be trained in reprocessing procedures of duodenoscopes. -Apply a robust microbiological surveillance system, reprocessing environments' cleanliness, and distal end of the duodenoscope design. -Develop schedules for routine inspection and periodic maintenance in accordance with the manufacturer's instructions. -Have a quality control program that includes sampling, microbiological culturing, and other monitoring methods for duodenoscopes. -Report all adverse events or complaints cases of using duodenoscopes to the SFDA. 	التوصيات
For Reporting	  	للإبلاغ