


## Safety Communication

### Counterfeit product

<b>Device/ Product Description:</b>	Vascular Prostheses
<b>Brand Name:</b>	UNI-Graft K DV
<b>Lot numbers/Serials:</b>	NA
<b>Manufacturer:</b>	Aesculap AG
<b>Problem:</b>	In the cause of an internal investigation, Aesculap AG became aware of product counterfeit of the UNI-Graft K DV vascular prostheses.
<b>Recommendation/Actions:</b>	<ul style="list-style-type: none"> <li>• Check all Uni-Graft K DV Vascular prostheses for correct packaging before use; prostheses that are not correctly packed and labelled as described below must not be used: <ul style="list-style-type: none"> <li>- Aesculap AG always delivers the products in double sterile packaging in cartons of 1 piece each. The counterfeit products are only single packed.</li> <li>- CE sign is missing on the counterfeit goods.</li> </ul> </li> <li>• You can read the complete Safety Alert that includes recommendations from ( <a href="#">HERE</a>).</li> </ul>
<b>Devices/Products photo:</b>	



You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative of your product for corrective action.

Healthcare Professionals should **report any adverse events** suspected to be associated with affected devices above (or other Medical Devices) to:

**National Center for Medical Devices Reporting.**

Medical Devices Sector

Saudi Food and Drug Authority

Postal Address: Saudi Arabia - Saudi Food and Drug Authority

4904 northern ring branch rd - Hitteen Dist.

RIYADH 13513 - 7148

Tel: +966 (11) 2038222 Ext: 2995, 2952

Fax: +966 (11) 2757245

Or

**Saudi Vigilance**

For latest published Recalls/Alerts, please visit ([NCMDR Website](#))

Sincerely,  
NCMDR Team

