

Safety Communication

رسالة سلامة

Foot and Ankle Instruments (Drill Tap and Countersink) tips may break

Device/ Product Description:	Foot and Ankle Inst. (Drill / Tap and Countersink)
Brand:	NORMED
Lot/Serial Number:	502015106 , 502015107 , 502015114 , 502015115 , 502015120 , 502015124 , 502015130 , 502015131 , 502015136 , 502015137 , 502015142 , 502015145 , 502015146 , 502015206 , ST502015206 , 502015207 , ST502015207 , 502015208 , ST502015208 , 502015211 , 502015212 , ST502015212 , 502015213 , 502015216 , ST502015216 , 502015217 , ST502015217 , 502015218 , 502015402 , ST502015402 , 502015403 , ST502015403 , 502015619 , ST502015619 , 502015620 , ST502015620 , ST502015621 , 502015623 , ST502015623 , 502015628 , ST502015628 , 502015629 , 502015630 , ST502015630 , 502015631 , ST502015631 , 502015634 , 502015635 , ST502015635 , 502015636 , 502015637 , 502015638 , 502015640 , ST502015640 , 502015650 , ST502015650 , 502015706 , ST502015706 , 503002041 , 503004177 , 503004341 , 503004342 , 503004351 , 503004352 , 503004353 , 503004541 , ST503004541 , 503004542 , ST503004542 , 28.66.110 , 28.66.111 , ST28.66.111 , 28.66.112 , 28.66.113 , ST28.66.113
Manufacturer:	Zimmer, INC
Problem:	Zimmer Biomet received a certain number of complaints reporting tip breakages. An investigation identifies that Normed Medizin Technik GmbH possibly manufactured certain lots from a different material than defined in the applicable specifications.


Recommendation/Actions:

1. Review this notice and ensure that affected personnel are aware of the contents.
2. If you have affected product at your facility, assist your Zimmer Biomet sales representative and quarantine all affected product. Your Zimmer Biomet sales representative will remove the affected product from your facility.

For more information, Please click [here](#).

Devices/Products photo:

NORMED Medizin-Technik GmbH: Packaging Information- **Affected by the Removal**



NORMED Medizin-Technik GmbH: Product Marking- **Affected by the Removal**



Normed Instruments- **Affected by the Removal**



Authorized Representative Details	AR name:	Medical Regulations Gate
	Assigned Contact Person:	Ameen almuslimani
	Mobile/Phone:	011 2934961
	Email:	ameen@mr-gate.com