

Safety Communication

رسالة سلامة

Detectors may fall - Discontinue system use until further notice

Device/ Product Description:	Gamma Camera Systems
AFFECTED PRODUCTS:	Forte, Forte JETStream, Forte JETStream upgrade, Forte JETStream AZ, Forte JETStream AZ upgrade, Diamond Select Forte, Diamond Select Forte JETStream
Manufacturer:	Philips Healthcare
Problem:	The manufacturer identified an issue affecting the mentioned products of gamma cameras that could result in either detector 1 or detector 2 falling unimpeded vertically to the end stops of its travel limit during detector radius movement. Possible contact with a patient could result in entrapment, serious injury, or death.
Recommendation/Actions:	<ul style="list-style-type: none">• Discontinue use of the system until further notice.• Inform those who need to be aware within your organization• Maintain this notice with your system Instructions for Use (IFU) until the correction is made to the system.• Communicate with the authorized representative to perform the required corrections.

Devices/Products photo:



**Authorized Representative
Details**

AR name:	Philips Healthcare Saudi Arabia Ltd.
Assigned Contact Person:	Mohammed AlSamhan
Mobile/Phone:	0558685708
Email:	sfda.sa.met@philips.com