

Safety Communication

رسالة سلامة

Philips Various Patient Monitoring Devices May Have Cybersecurity Vulnerabilities

| | | |
|-------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|
| Device/ Product Description: | Central Monitoring Systems, Patient Monitors | |
| Manufacturer: | Philips | |
| Problem Summary: | <p>Cybersecurity vulnerabilities that could allow an attacker with a low skill level unauthorized access, potentially leading to interrupted monitoring and collection of access information and/or patient data.</p> <p>For more information, please check Here.</p> <p>If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through: NCMDR Vigilance system 19999 unified call center</p> | |
| | AR name: | Philips Healthcare Saudi Arabia Ltd. |
| | Assigned Contact Person: | Alawi Alkhadrawi |
| | Mobile/Phone: | 0558685708 |
| | Email: | sfda.sa.met@philips.com |