

Medical Devices Sector
National Center for Medical Devices Reporting

رسالة سلامة
Safety Communication

To: Healthcare Providers	إلى: مقدمي الرعاية الصحية
---------------------------------	---------------------------

Title:	Risk of catheter breakage	العنوان:
Medical Device Description	- IV Cannula G: 20*32mm non port - IV Cannula G: 22*25mm non port	اسم ووصف الجهاز/المستلزم الطبي:
Manufacturer:	Saudi Mais Co, for Medical Products	اسم المصنع:

Potential /Associated Risks	There was insertion of IVC in Jugular vein and also re-insertion of needle in catheter which lead to breakage of catheter.	المخاطر المحتملة المرتبطة بالجهاز أو المستلزم الطبي
Recommendations	The IV catheter is always recommended to use in Peripheral sites and avoid re insertion of Needle, because it may cause breakage of catheter	التوصيات

For Reporting:	  	لإبلاغ:
-----------------------	--	---------